




ADMINISTER  
MODULAR  
NUTRITION  
AT THE SPEED  
OF NURSE.



# 1. OPEN



# 2. ATTACH



# 3. SQUEEZE



Care in the ICU is fast paced and often requires multiple lines of intervention. Nutrition support plays an important role in patient care but can be labor intensive and time consuming. ENFit®-compatible modular nutrition allows for easy intervention that is efficient and effective.

 More than nutrition,  
**Medtrition®**

## ProSourceTF20- Meeting protein goals has never been easier.

 **60.5%**

Studies show ICU patients receive only 60.5% of their protein requirements over 4 days and 66.7% over a 12-day period.\* ASPEN guidelines recommend 1.2g-2.0g per 1 kg body weight for ICU patients.


Protein is the most important macronutrient in the ICU, yet many patients don't get enough. ProSource TF20 makes administering protein easier than ever.

- Nurses can administer quickly and easily
- Water thin
- Direct administration, no mixing or stirring
- Good GI tolerance
- Each 60 ml sachet contains 20 grams of protein
- Helps meet ASPEN critical care guidelines for protein

### ADMINISTER PROTEIN MORE EFFICIENTLY.

Most enteral formulas don't provide the required amount of protein to meet the needs of critically ill patients. Even costly specialty formulas fall short. ProSource TF20 helps meet protein goals with fewer packets, saving money and time. With these products nurses also will use fewer cups and syringes.

\* Nicolo M, Heyland DK, Chittams J, Sammarco T, Compher C. Clinical Outcomes Related to Protein Delivery in a Critically Ill Population: A Multicenter, Multinational Observation Study. JPEN J Parenter Enteral Nutr. 2016 Jan;40(1):45-51. doi: 10.1177/0148607115583675. Epub 2015 Apr 21. PMID: 25900319.



**100**  **PDCAAS Score**

ProSource TF 20 rates the highest PDCAAS score. PDCAAS evaluates a food's protein quality by comparing its amino acid composition to what our bodies can use.

## Banatrol TF— The next evolution in managing diarrhea

Banatrol TF solidifies stool without causing constipation. For the first time, nurses can easily administer an all-natural solution to relieve the burden of cleaning diarrhea and spare patients the emotional and physical harm it causes.

- Nurses can administer quickly and easily
- Liquid fiber formula
- Direct administration, no mixing or stirring
- Made with banana flakes, soluble fiber and prebiotics
- Promotes long-term gut health
- Helps meet ASPEN critical care guidelines for diarrhea intervention

### BANATROL TF—MANAGE DIARRHEA FROM DAY ONE.

Banatrol TF can be administered when a patient has diarrhea and prophylactically to prevent diarrhea as soon as a patient begins tube feeding. Banatrol is used by hospitals and patients across the country.

### HELPS MEET ASPEN GUIDELINES.

The Banatrol TF formula works with a patient's digestive system. It helps meet 2016 ASPEN critical care guidelines for diarrhea intervention – 10-20 grams of fermentable soluble fiber to be given in divided doses over 24 hours for those who have signs of diarrhea.\*

### HOW IT WORKS

Banatrol TF contains resistant starch to feed the good bacteria in the gut and pectin to help absorb water and improve the stool consistency. Clinical studies show the benefits of soluble fiber in nutrient absorption, gut motility and stool characteristics. Prebiotic fiber plays an important role in the microbial balance of the GI system.\*

\* McClave, S.A., Taylor, B.E., Martindale, R.G., Warren, M.M., Johnson, D.R., Braunschweig, C., McCarthy, M.S., Davanos, E., Rice, T.W., Cresci, G.A., Gervasio, J.M., Sacks, G.S., Roberts, P.R., Compher, C., and (2016), Guidelines for the Provision and Assessment of Nutrition Support Therapy in the Adult Critically Ill Patient. Journal of Parenteral and Enteral Nutrition, 40: 159-211. <https://doi.org/10.1177/0148607115621863>.




 **More than nutrition, Medtrition®**




# READY FOR THE NEW STANDARD

ENFit® is a global change to make all tube-feeding devices specific to tube feeding. ProSource TF20 and Banatrol TF are ENFit®-enabled to make administering easy whether your facility is ENFit®-enabled or not.




For traditional legacy tubes, nurses need only


1) Open



2) Squeeze package into a cup



3) Syringe into tube.



## Nutrition isn't one size fits all

Standard enteral formulas often fall short in protein provision especially when challenged with the higher-end needs of some patients. Medtrition products work in unison with standard enteral formulas so dietitians can customize and nurses can administer care specific to the patient's needs.

## MORE THAN NUTRITION

Medtrition develops groundbreaking medical foods used in nursing homes and hospitals across the U.S. and in over 30 countries. Our science and products enable dietitians, nurses and facilities to stay on the leading edge of healing through nutrition.

To learn more, contact your sales representative.  
Visit [Medtrition.com/TF20](http://Medtrition.com/TF20) or [Medtrition.com/BanatrolTF](http://Medtrition.com/BanatrolTF)  
Or call 877-271-3570 for product or company details.

More than nutrition,  
**Medtrition®**